

Easy and Fast



NeoPlex™

Innovative Multiplex  
Real-time PCR Assay

NeoPlex™ SARS-CoV-2

# COVID-19 Detection

US EUA Approved  
CE-IVD Marked



## Single-Tube Real-time RT-PCR Assay



## Dual Targets recommended by US CDC and WHO

- N and RdRp genes



## Fast Turnaround Time

- Within 2 hours analysis after nucleic acid extraction



## Ease of Use


- Automated analysis with NeoPlex™ Viewer software
- One Step RT-PCR without separate cDNA synthesis step
- Whole process monitoring with internal control



## Compatible Real-Time PCR instruments

- CFX96™ Real-time PCR (Bio-Rad)
- Applied Biosystems™ 7500 (Fast) Real-time PCR (Thermo Fisher)
- Applied Biosystems™ 7500 Quantstudio5 (Thermo Fisher)
- Gentier 96 Real-time PCR (TianLong)
- Rotor-Gene Q 5plex HRM (Qiagen)
- LightCycler 480 (Roche)

## NeoPlex™ Viewer



## Specimens

Nasopharyngeal or Oropharyngeal swab, Sputum, Saliva, Bronchoalveolar lavage fluid(BAL)

## Series Products

NeoPlex™ COVID-19 FAST : RdRp, N

NeoPlex™ COVID-19 Variant : RdRp, Spike H69-V70 del., Spike L242-L244 del., ORF1a S3675-F3677 del.

NeoPlex™ FluCOVID : SARS-CoV-2 (RdRp, N), Flu A, Flu B

## Ordering Information

Cat No.	Product	Description	Method	Size
NR05A	NeoPlex™ COVID-19	Simultaneous Detection of N and RdRp genes for SARS-CoV-2 in a Single Test	Multiplex Real-time RT-PCR	96 Tests/Kit
NR12A	NeoPlex™ COVID-19 FAST	Fast and Simultaneous Detection of N and RdRp genes for SARS-CoV-2 in a Single test	Multiplex Real-time RT-PCR	96 Tests/Kit
NR11A	NeoPlex™ COVID-19 Variant	Simultaneous Detection of SARS-CoV-2 and its Variants in a Single Test	Multiplex Real-time RT-PCR	96 Tests/Kit
NR07A	NeoPlex™ FluCOVID	Simultaneous Detection of SARS-CoV-2 and Influenza A/B in a Single Test	Multiplex Real-time RT-PCR	96 Tests/Kit

### GeneMatrix Inc.

KoreaBio Park B Bldg. 8F  
 700 Daewangpangyo-ro, Bundang-gu  
 Seongnam-si, Gyeonggi-do, 13488 South Korea  
 Tel : +82-31-628-2020  
 Email : sales@genematrix.net